



METECH
— INTERNATIONAL —

Professional Eye
Instruments Range



Ask For it!

info@metech-international.com

Our Products Are Guaranteed And Meet The Relevant
Worldwide Standards And Our Products Are
Manufactured Entirely Under In Factory.

COMPANY PROFILE



Metech International Based In Sialkot And Oldest Manufacturer And Exporter Of Vast Rang Of Instruments. It comprises of an integrated manufacturing facility, employing skilled craftsmen to produce broad range of professional Instruments. Its manufacturing process is running under the supervision of expert technicians.

MANAGEMENT

Skilled employees are its main source for innovation and development of new products. They are working in safe and open environment where every worker contributes continuously to expand the company business. Qualified persons are hired to make use of the progressive technology. The company is also providing active support to its employees to improve their skills, creativity power and social development.

Our Aim

Metech International & Company has earned the credit to get the entire satisfaction of the Customers. By providing

Quality Products
Fair Prices
Timely Deliveries
Product Range

Metech Intl produces different types of Instruments. Our Main line of the business is General Surgical Instruments, Dental, Plastic Surgery Instruments, ENT Instruments, Gynecology & Ophthalmic Instruments. Each Instrument manufactured by **Metech Intl** has personally been supervised and monitored to make sure that the Instruments are perfect. In other words, continuous research and development has allowed this company to make superior quality products.


Quality Standards

Metech Intl has achieved the certification of Sialkot Chamber of Commerce, Surgical Instruments Manufacturer Association Of Pakistan, NTN, CE Mark, ISO 13485-2016. These quality assurance standards are maintained throughout the different production process by its modern testing devices and trained staff which is the main hallmark to maintain the quality of its products.



Ask For it!

info@metech-international.com
Our Products Are Guaranteed And Meet The Relevant
Worldwide Standards And Our Products Are
Manufactured Entirely Under In Factory.



CORPORATE GOAL

To be the market leader in manufacturing world's finest Surgical & Dental Instruments, at the most competitive prices.

Never compromise on Quality!

During the initial years of our formation, we worked hard to precisely understand the requirements of our respected customers, which has gradually improved over time. As our guiding philosophy, we put the concepts of "Customer First" and "Quality First" into practice, to respond to the expectation of our customers.

QUALITY ASSURANCE

Our commitment to the quality ethics and performance of all our instruments is demonstrated by our accreditation to the International quality assurance standards of ISO 9001:2008, ISO 13485:2003, cGMP & CE.

GUARANTEE

Our guarantee is your complete satisfaction.

We take great pride in our reputation for manufacturing "consistent" quality instruments as per our customer's requirements. We offer our customers a complete guarantee when it comes to the workmanship and materials of all our products.

R&D

Our company take R&D as a continuing driving force for the research and development of our Instruments. Our company is well equipped with in-house R&D department, with state of the art machinery, where every piece is put through stringent tests before its final assembling. It provides us ample scope to adapt the needs of our valued clients.

CHILD LABOR POLICY


As good corporate citizens and most importantly as a matter of social responsibility, we are committed to the ideas of protecting children from the exploitations of child labor. We strongly believe that the future growth of children, of their communities and the country as a whole is best served through education, not child labor.



REPAIR & RESTORATION

We offers complete surgical and dental instruments repair and restoration service, utilized in all kinds of operating procedures at very competitive rates.

Our services include, but are not limited to, complete set restoration, instrument customization, gold plating, laser/chemical etching, coding, final finishing, Re-Polishing and most importantly the functional performance of an instrument and its aesthetic appearance.




OEM MANUFACTURING

We not only manufactures fine quality surgical instruments, but we also provide complete modification and customization of instrument design according to its specific usage. All OEM instruments are carefully hand crafted with details to guarantee quality and exact specifications.

We can work pretty much on anything, from physical samples to hand drawn sketches, high precision 3D CAD models and drawings.

Our experience and knowledge in hand crafting surgical and dental instruments, combined with modern engineering capabilities, allow us to adopt the finest techniques to apply to each product, to accomplish the outcomes our customers require.



FORCEPS	06
MARKERS CLIPPERS	17
SCISSORS	18
FORCEPS	22
SPATULAS	23
NUCLEUS MANIPULATORS	24
SPATULAS LOOPS	25
PUNCHES POLISHERS	26
SPATULAS MANIPULATORS	27
LOOPS	28
SPOONS IRIS RETRACTORS	29
NUCLEUS DIVIDERS SPOONS	30
IRIS HOOKS	31
IOL HOOK	32
HOOKS	33
NEEDLE HOLDERS	35
SPECULUMS	39
RETRACTORS	45
KNIVES	47
KNIVES SPUDS	48
CANNULAS	49
BLADE BREAKERS	50
I-A CANNULAS	51
CYSTOTOMES	52
BLADE HANDELS	53
CAPSULE POLISHERS	54
DISPOSABLE INSTRUMENTS	56
CANNULAS	57
CATRARACT SET	59
INSTRUMENTS SET	60
STERILIZING CASE	85



Gill
 Iris Forceps, curved
 criss-cross serrated tips
 1x2 teeth
FS-1001



Botvin
 Iris Forceps
 for peripheral iridectomies
 serrated tips
 1x2 teeth
FS-1002



Gills-Colibri
 Forcep
 very fine pointed tips
FS-1003



Colibri
 Forcep, very delicate
 1x2 teeth, 0.12mm
FS-1004



Harms-Colibri
 Forcep, very delicate
 1x2 teeth, 0.12mm
 with tying platform
FS-1005



Dyson-Colibri
 Forcep, very delicate
 1x2 teeth, 0.12mm
 extremely thin shanks
 short tips with tying platform
FS-1006



Colibri

Forceps, standard
 1x2 teeth, 0.4mm
 with tying platform
FS-1007



Polack

Corneal suturing forceps
 1mm spread, very delicate
 1x2 teeth, 0.12mm
FS-1008



Hofmann-polack

Corneal suturing forceps
 with guide pin, 1mm spread
 1x1 pierse teeth, 0.2mm

FS-1009



Pierse-colibri

Corneal forceps
 very delicate
 1x1 pierse teeth, 0.1mm

FS-1010



Pierse

Corneal forceps
 very delicate
 1x1 pierse teeth, 0.1mm
 straight
 curved

FS-1011



Castroviejo-Colibri

Forceps very delicate
 1x2 teeth, 0.12mm

FS-1012





Girards

Corneoscleral Forceps
 1x2 teeth, 0.12mm
 with tying platform
 right handed
 left handed

FS-1013



Maumenee-Colibri

Corneoscleral Forceps
 very delicate
 1x2 teeth, 0.12mm

FS-1014



Tennant-Colibri

Corneoscleral Forceps
 very delicate
 1x2 teeth, 0.12mm

FS-1015



Bores

Corneoscleral Fixation Forceps
 U-shaped, 3mm spread
 1x2 teeth at each tip,
 0.12mm

FS-1016



Rubman Bore

Corneoscleral Fixation Forceps
 U-shaped, 3mm spread
 1x2 teeth at each tip,
 0.12mm

FS-1017



Arrowsmith

Corneoscleral Fixation Forceps
 with lock
 U-shaped, 8mm spread
 1x2 teeth at each tip,
 0.12mm

FS-1018



Dressing

Forceps, serrated
delicate
straight
curve

FS-1019



Tissue

Tissue Forceps
1x2 teeth
delicate
standard

FS-1011



Tissue

Forceps, 1x2 teeth
delicate
straight
curved

FS-1021



Lester

Fixation Forceps
straight
1x2teeth
2x3teeth

FS-1022



Barken

Fixation Forceps
1x2 teeth

FS-1023





Katena

Corneal Fixation Forceps
 with lock
 U-shaped, 13mm spread
 1x2 teeth at each tip,
 0.12mm
 straight
 angled
FS-1024



Bishop-Harmon

Tissue Forceps
 1x2 teeth
 delicate
 standard
FS-1025



Bishop-Harmon

Dressing Forceps
 criss-cross serrate
FS-1026



Sauer

Fixation Forceps
 1x2 angled teeth
FS-1027



Manhattan Eye & Ear

Fixation Forceps
 1x2 angled teeth
 with tying platf
FS-1028



Farkas-Bracken

Fixation Forceps
1x2teeth, delicate

FS-1029



Kuhnt

Fixation Forceps
1x2 curved teeth

FS-1030



Elschnig

Fixation Forceps
1x2 angled teeth

FS-1031



Nugent

Superior Rectus Forceps
angled, 1x1 teeth

FS-1032



Troutman

Superior Rectus Forceps
angled, 1x2 teeth

FS-1033



Tennant-Troutman

Superior Rectus Forceps
curved shanks, 1x2 teeth


FS-1034



Graefe

Fixation Forceps
4.5mm wide jaws
without lock
with lock

FS-1035






Tennant

Tying Forceps
 with guide pin
 extra delicate
 smooth jaws with
 6mm long platform
 for 9-0 to 11-0 sutures
 straight
 curved

FS-1036



Harms

Tying Forceps
 smooth jaws
 for 7-0 to 9-0 sutures
 straight
 curved

FS-1037



Fine

Tying Forceps
 smooth curved jaws
 for 7-0 to 9-0 sutures

FS-1038



Barraquer

Cilia Forceps
 lightly angled
 smooth curved jaws

FS-1039



Katena

Cilia Forceps
 slanted jaws
 gold plated

FS-1040



Gradle

Cilia Forceps
 smooth, squared off
 positive gripping jaws

FS-1041



Lambert
 Chalazion Forceps
 round, 8mm diameter
FS-1042



Hunt
 Chalazion Forceps
 round, 8mm diameter
FS-1043



Baird
 Chalazion Forceps
 oval, small
FS-1044



Francis
 Chalazion Forceps
 oval, medium
FS-1045




Desmarres
 Chalazion Forceps
 oval, large
FS-1046



Wies
 Chalazion Forceps
 grooved upper jaw
FS-1047



Lordan
 Chalazion Forceps
 triangular jaws
FS-1048





Jewelers
 Forceps, delicate
FS-1049



Jewelers
 Forceps, delicate
FS-1050



Jewelers
 Forceps, extra delicate
FS-1051



Osher
 Iris Forceps
 very fine textured
 jaws
FS-1052



Jewelers
 Forceps, extra delicate angled
FS-1053



Jewelers
 Forceps, curved
FS-1054



Suture Removal Forceps
 Very fine pointed tips
 with carbide impregnated jaws
FS-1055



Serrefine
 Small, straight
FS-1056



Hartman
 Mosquito Forceps
 extra delicate
 3 3/4" (95mm) long
 straight
 curved
FS-1057



Halsted
 Mosquito Forceps
 delicate
 5 3/4" (145mm) long
 straight
 curved
FS-1058



Backhaus
 Towel Clamp
 3 1/2" (90mm) long
 5 1/4" (135mm) long
FS-1059



Baby Jones
 Towel Clamp
 cross-action
 2 1/4" (55mm) long
FS-1060





Shepard II
 Lens Holding Forceps
 with nudge notch
 for A/C IoLs
 with lock
 without lock
FS-1061



Jameson
 Muscle Forceps
 with slide lock
 4teeth, child size
 right
 left
FS-1062



Jameson
 Muscle Forceps
 with slide lock
 6teeth, adult size
 right
 left
FS-1063



Berke
 Ptosis Forceps
 with slide lock
 longitudinally grooved jaws
 20mm long jaws
 27mm long jaws
FS-1064




Erhardt
 Lid Forceps
FS-1065





Anis

Radial Marker
 40mm center opening
 12.5mm O.D.
 3 blades 10 blades
 4 blades 12 blades
 6 blades 16 blades
 8 blades



MC-1071

Ains

Suture Placement Marker
 for corneal grafts, 8 radial blades
 with 8mm center opening



MC-1072

Deitz

Incision Depth Gauge
 for corneal incisions
 0.5mm
 0.55mm



MC-1073

Thornton

Limbal Incision Ruler
 with marking lines
 lines spaced 1mm apart*



MC-1074

Lindstrom


Small Incision Marker
 for cataract incision widths of
 3, 4, and 5mm at 1, 2, or 3mm
 away from the limbus*



MC-1075

Castroviejo

Caliper, 20mm spread
 straight
 curved



MC-1076

Osher

Incision Caliper
 for internal measurement of
 cataract incisions from 1 to 10mm
MC-1077



Castroviejo
 Keratoplasty Scissors
 angled to side
 blunt tips
 small blades
SC-1066



Castroviejo
 Keratoplasty Scissors
 angled to side
 blunt tips
 medium blades
SC-1067



McPherson-Westcott
 Conjunctival Scissors
 curved, blunt tips
 small blades
SC-1068



Westcott
 Tenotomy Scissors
 curved, blunt tips
 medium blades
SC-1069



Westcott
 Tenotomy Scissors
 curved, blunt tips
 standards blades
 right
 left
SC-1070



Castroviejo
 Corneoscleral Scissors
 small blades
 right
 left
SC-1078



Castroviejo
 Corneoscleral Scissors
 Medium blades
 Right
 Left
SC-1079



Troutman-Castroviejo
 Corneoscleral Scissors
 with stop, small blades
 Right
 Left
SC-1080





Vannas
 Capsulotomy Scissors
 sharp tips
 angled to side
SC-1081



Gills-Vannas
 Capsulotomy Scissors
 extra thin long blades
 sharp tips
 medium blades
SC-1082



Gills-Vannas
 Capsulotomy Scissors
 extra thin long blades
 sharp tips
SC-1083



Gills-Vannas
 Capsulotomy Scissors
 extra thin long blades
 sharp tips
SC-1084



McPherson-Vannas
 Iris Scissors
 sharp tips
 curved
SC-1085



Katena-Vannas
 Scissors, extra thin
 7mm long blades
 sharp points
 straight
SC-1086



Eye Scissors
 Standard
 pointed tips
 straight
 curved
SC-1093



Eye Scissors
 Standard, large rings
 pointed tips
 straight
 curved
SC-1094



Stitch Scissors
 Ribbon type
 needle points
 light curved
SC-1095



Shaepard

Lens holding Forceps
 with serrated lower jaw
 gently curved shanks
 without lock

with lock

FC-1087



Sinskey

Lens holding Forceps
 delicate T-shaped lower
 jaw for 3-points fixation

FC-1088



Kraff

Lens holding Forceps
 delicate curved jaws
 to grasp edge of lens optic

FC-1089



Anis

Lens Inserting Forceps
 very delicate curved tips
 to grasp IOL haptics.
 Angled shanks for
 comfortable hand position

FC-1090



Clayman

Lens holding Forceps
 delicate angled jaws
 gently curved shanks
 without lock


with lock
FC-1091



Phaco-Folder II

For folding and inserting
 soft IOLs through a
 small incision

FC-1092



Barraquer
 Corneoscleral brush
 Sable with plastic handle
SP-1096



McPherson
 Iris Spatula
 1mm wide
 malleable stainless steel
SP-1097



Culler
 Iris spatula
 Sterling silver
 Blade
 1mm wide
 2mm wide
SP-1098



Wheeler
 Double-Ended Spatula
 Lightly curved
SP-1099



Barraquer
 Iris Spatula
 Round, smooth
 0.5mm dia.
 0.25 mm dia.
SP-1100



Gimbel
 Nucleus Spatula
 Angled, round, smooth
 0.25mm dia.
SP-1101



Knolle
 Lens Nucleus Spatula
 Malleable stainless steel
 Flat, 0.5mm wide
SP-1102





Bechert
 Lens Nucleus Rotator
 Blunt forked tip
 Angled at 7mm
NM-1103



Jaffe-Bechert
 Lens nucleus Rotator
 Blunt forked tip
 Angled at 9mm
NM-1104



Clayman
 Nucleus Rotator
 0.5mm wide with
 V-notched tip
NM-1105



Drysdale
 Nucleus Manipulator
 Paddle shaped tip
NM-1106



Pisacano
 Nucleus Rotator
 Flat blunt tip
NM-1107



Siepser
 Endocapsular Controller
 Trumpet shaped tip
NM-1108



Tennant
 Nucleus Manipulator
 Ball shaped tip
NM-1109



Elsching
 Cyclodialysis Spatula
 Flat, 1mm wide blade
SL-1110



Castroviejo
 Synechiaie Spatula
 semi-sharp tips
 0.5mm wide
 double-end
SL-1111



Castroviejo
 Cyclodialysis Spatula
 double-end
 0.25mm wide
 0.5mm wide
SL-1112



Kimura
 Platinum Spatula
SL-1113



Lewis
 Lens Loop
 Small
 Large
SL-1114



New Orleans
 Lens Loop
SL-1115



Wilder
 Lens Loop
 serrated
SL-1116



Gass

Corneoscleral punch
 1.5mm dia. Tip
 rotatable for
 Convenient positioning
PP-1117



Kelly

Descemet's
 Membrane Punch
 1.0mm dia. tip
PP-1118



Meyerhoefer

Chalazion Curette
 Size 0, 1.5mm dia.
 Size 1, 1.75mm dia.
 Size 2, 2.25mm dia.
 Size 3, 2.5mm dia.
 Size 4, 3.5mm dia.

PP-1119



Gills

Capsule Polishing Curette
 1.75mm dia. Cup
 with sandblasted base
 gentel backward angle
PP-1111



Nightingale

Capsule Polishing Curette
 1.75mm dia. Open ring
 with semi-sharp edge
PP-1121



Anis

Capsule Scrubber
 Gently curved flexible shank
 With roughened tip
 Ball tip
 Disc tip
 Reverse ball tip

PP-1122



Lippman

Sideport Nucleus Spatula
 0.5mm dia. Round
 blunt forked tip
SM-1123



Girard

Synechiaie Spatula
 Semi-sharp tip
 0.5mm wide
SM-1124



Obstbaum

Lens Spatula
 Angled blade
 0.25mm wide
SM-1125



Hirschman

Lens Spatula
 Flat blade, 0.5mm wide
 With small notches for
 Manipulating IOLs
SM-1126



Jaffe

Lens Spatula
 Flat, smooth
 0.5mm wide
SM-1127



Lester

IOL Manipulator
 For Lenses with or
 Without hole

Straight, with guard
 angled

SM-1128





Kirby
 Lens Loop
 And Muscle Hook
L-1129



Gills-Welsh
 Lens Loop and Probe
 with serrations
L-1130



Morrison
 Lens Loop
 And Probe
 Very thin
L-1131



Simcoe
 Lens Nucleus Loop
 Fine serrations, strong curve
 Double-end
L-1132




Kansas
 Lens Loop, 3mm wide
 Serrated, for small incision
 Nucleus fragment
 Removal
L-1133




Keener-Arit
 Lens Loop
 3mm wide, serrated
 strong curve for
 small incision nucleus
 fragment removal
L-1134



Keener
 Nucleus Divider
 With two pre-formed
 Snare wire loops
 Complete
 Wire loops(pk/3)
L-1135



Bunge
 Evisceration Spoon
 Small
 Large
SI-1136



Wells
 Enucleation Spoon
SI-1137




Rizzuti
 Iris Retractor
 And Lens
 Expressor
 2mm wide lip
SI-1138



Azar
 Iris Retractor
 Autoclavable plastic
 Left
 Right
SI-1139



Graether
 Collar Button
 Micro iris retractor
 Straight, with guard
 angled
SI-1140



Kansas

Nucleus Trisector
 For manual small incision surgery
 For right handed
 Surgeons
 For left handed
 Surgeons
NS-1141



Kansas

Nucleus Bisector
 For manual small incision surgery
NS-1142



Kansas

Nucleus Vectis, solid
 Used as a "cutting block"
 For Kansas nucleus
 Bisector and trisector
 For right handed surgeons
 For left handed surgeons
NS-1143



Culler

Lens Spoon
 Sterling Silver
NS-1144



Barraquer

Lenticular Spoon
 For corneal button transfer
NS-1145



Paton

Spatula and Spoon
 For corneal button transfer
NS-1146




Graether
 Collar Button
 Micro iris retractor
 And posterior Capsule polisher
 With side irrigation
 Straight
 Angled
IH-1147



Bonn
 Micro Iris Hook
 Blunt, round handle
 With protective Guard
IH-1148



Bonn
 Micro Iris Hook
 Blunt, flat handle
IH-1149



Jaffe-Knolle
 Iris Hook, very blunt
 Small, angled
IH-1150



Knolle
 Iris Hook, very blunt
 Small, straight
IH-1151




Shepard
 Iris Hook
 Small, blunt
IH-1152






Sinskey

Lens Manipulating Hook
 With round handle
 0.25mm dia. Blunt tip
 straight, with guard
 angled
IH-1153




Sinskey II

Lens manipulating Hook
 With round handle
 0.2mm dia. Blunt tip
 straight
 angled
IH-1154



Sinskey II

Lens Manipulating Hook
 With round handle
 0.2mm dia. Blunt tip
 straight, with guard
 angled
IH-1155



Lewcky

Lens manipulating Hook
 0.15mm dia. Blunt tip
 angled
IH-1156



Hirschman

Iris Hook
 Blunt, angled
IH-1157



Maltzman-Fenzl

Lens Manipulating Hook
 0.15mm dia. V-shaped tip
 round handle
 straight, with guard
 angled
H-1158



Galand

Lens Manipulating Hook
 angled
H-1159



Jaffe-Maltzman

Lens Manipulating Hook
 0.15mm dia. V-shaped tip
 angled
H-1160



Osher

Y-Hook
 for manipulating IOLs
 straight, with guard
H-1161



Fixation Forceps

double, sharp
 small
 large
H-1162




Twist Hook

for scleral fixation
 right
 left
H-1163



Anis

Scleral Twist Fixation Cone
 2.5mm dia. With three
 teeth, with guard
H-1164





Kuglen

Iris Hook
And Lens Manipulator
Push-pull model
Straight
Angled

H-1165



Kuglen

Iris Hook
And Lens Manipulator
Clover leaf model
Straight

H-1166



Kuglen

Iris Hook
And Lens manipulator
Push-pull model
With round handle
Straight, with guard
Angled

H-1167



Rentsch"Boat Hook"

Lens Manipulator
This multi-purpose
Instrument may be used as
An iris/capsule retractor or
Pusher/manipulator for IOLs
With and without holes
Straight
With guard

H-1168




Fenzl

Lens Manipulating Hook
0.15mm dia. T-shaped tip
straight,
with guard
angled

H-1169



McPherson
 Needle Holder
 very delicate jaws
 straight, with lock
 straight, without lock
 curved, with lock
 curved, without lock
NH-1170



Jaffe
 Needle Holder
 delicate conical jaws
 straight, with lock
 curved, with lock
 curved, without lock
NH-1171



Castroviejo
 Needle Holder
 delicate jaws
 straight, with lock
 straight, without lock
 curved, with lock
 curved, without lock
NH-1172





Castroviejo
 Needle Holder
 standard jaws, with lock
 straight
 curved
NH-1173



Kalt
 Needle Holder
 straight smooth jaws
 delicate
NH-1174



Drews
 Needle Holder
 extra delicate jaws
 straight, with lock
 straight, without lock
 curved, with lock
 curved, without lock
NH-1175



Barraquer
 Needle Holder
 short model
 very delicate jaws
 straight, with lock
 straight, without lock
 curved, with lock
 curved, without lock
NH-1176



Anis

Needle Holder
 extra delicate jaws
 straight, with lock
 straight, without lock
 curved, with lock
 curved, without lock
NH-1177




Katena

Needle Holder
 extra delicate jaws
 straight, with lock
 straight, without lock
 curved, with lock
 curved, without lock
NH-1178



McInyre

Fish Hook Needle Holder
 extra delicate jaws
 straight, with lock
NH-1179



Troutman

Needle Holder
 delicate, strongly curved jaws,
 with lock
 without lock
NH-1180



Barraquer

Needle Holder
 very delicate curved jaws
 with lock
 without lock
NH-1181



Barraquer

Needle Holder
 delicate jaws
 straight, with lock
 straight, without lock
 curved, with lock
 curved, without lock
NH-1182





Barraquer

Needle Holder
 short model
 standard jaws, curved
 with lock
 without lock
NH-1183




Barraquer
 Wire Speculum
 small
 large
 large, heavy wire
SP-1184




Barraquer
 Wire Speculum
 with solid blades
 small
 large
SP-1185





Grandon-Barraquer
 Wire Speculum
 with suture notches or holes
 to obtain additional traction
 for the upper lid
 wire blades
 solid blades
SP-1186




Feaster-Barraquer
 Wire Speculum
 wide, solid blades
SP-1187



McIntyer
 Wire Speculum with
 gently curved spring.
 V-shaped blades help to retain
 disposable drape under lids
 close wire
 open wire
SP-1188



Kartz-Barraquer
 Wire Speculum
 with open blades
 for unobstruted access
 to the cataract incision
 with phaco or I/A tip
 heavy wire, large
 standard wire, large
 standard wire, small
SP-1189



Kartz-Barraquer
 Wire Speculum with wings
 for retracting surgical
 drape while allowing
 unobstruted access
 to the cataract incision
 with phaco or I/A tip
 standard wire, small
SP-1190



Tennant
 Lid Retractor
 may be clamped to surgical
 drape with suture or elastic band
 pair
SP-1191



Sauer
 Lid Speculum
 for infants
 for premature infants
SP-1192



Infant
 wire Speculum
 with gently curved springs
 4mm wide blades
SP-1193



Cook
 Eye Speculum with locking screw
 solid, small #1
 solid, medium #2
 solid, large #3
 fenestrated, medium #2
 fenestrated, large #3
SP-1194



Castroviejo

Eye Speculum
 with finely controlled opening
 and closing mechanism
 small #1
SP-1195



Castroviejo

Eye Speculum
 with finely controlled opening
 and closing mechanism
 medium #2
SP-1196




Fishkind-castroviejo

Lid Speculum with open blades
 for use with phaco I/A tip
SP-1197



McPherson

Eye Speculum
 solid blades with
 vertical suture posts
SP-1198



Williams
 Eye Speculum
 open blade, small
SP-1199



Williams
 Eye Speculum
 open blade, large
SP-1200



Lancaster
 Eye Speculum
 solid blade
SP-1201



Lester-Burcs
 Eye Speculum
 solid blade with
 flanges for retaining
 eye lashes
SP-1202



Meaumennee-park

Eye Speculum
 with canthus hook,
 blade locking screws may
 be positioned vertically
 or horizontally
 solid blades
 fenestrated blades

SP-1203

Kuglen


Eye Speculum
 with canthus hook,
 solid blades with
 swivel ball joints
 and locking screws

SP-1204

Guyton-park

Eye Speculum
 with vertical suture posts
 and canthus hook
 solid blades
 fenestrated blades

SP-1205



SPECULUMS


Smirmaul

Eye Speculum
individually adjustable
blade for elevating lids
solid blades
wire blades
SP-1206




Schott

Eye Speculum
with rotatable adjustable
for precise control of
blades elevation
wire blades
SP-1207



Flieringe

Scleral Fixation Ring
set of 8 sizes(15-22mm)
14mm dia.
15mm dia.
16mm dia.
17mm dia.
18mm dia.
19mm dia.
20mm dia.
21mm dia.
22mm dia



SP-1208



Jaeger

Lid plate
 stainless steel
RE-1209



Wright

Fascia Needle
 with 1mm x 6mm
 oval hole
 14cm long
RE-1210



Knapp

Lacrimal Retractor
 4-prong
 blunt
 sharp
RE-1211



Stevenson

Lacrimal Sac Retractor
 3 x 3 blunt prongs
 self-retaining
 adjustable spread
RE-1212




Stevenson

Lacrimal Sac Retractor
 solid blades with
 serrated edge
 self-retaining
 adjustable spread
RE-1213



Agricola

Lacrimal Sac Retractor
 3 x 3 sharp prongs
RE-1214



Schepens
 Forked Orbital Retractor
RE-1215




Arruga
 Globe Retractor
RE-1216



Jaffe
 Wire Lid Retractor
 may be clamped to
 surgical drape with
 suture or elastic band
 pair, large
 pair, small
RE-1217



Desmarres
 Lid Retractor
 solid blade
 size 0
 size 1
 size 2
 size 3
RE-1218



Jaeger
 Lid Plate
 clear plastic
RE-1219



Ziegler

Knife Needle
 3mm
 5mm
KN-1220



Girard-swan

Knife Needle
 1.5mm cutting edge
 5mm
KN-1221



Wheeler

Discission Knife
 18mm cutting edge
KN-1222



Emery

cystotome
 double-edge
 for push and pull cutting
 blade, straight
 blade, angled
 chuck handle
KN-1223




Stiletto Blade

For micro-vitro-retinal
 stab incisions, lancet
 shaped, 19 gauge
 1.2mm wide
 blade
 chuck handle
KN-1224



Paufique

Graft Knife
 angled tip
 curved cutting edge
KN-1225





Tooke

Coeneal Knife
 straight blade with
 curved cutting edge
KS-1226



Gill

Coeneal Knife
 curved cutting edge
 straight
 curved
KS-1227



Golf-Club

Foreign Body Spud
 angled, lightly rounded tip
KS-1228



Francis

Foreign Body Spud
 angled, flag shaped tip
KS-1229



Ellis

Foreign Body Spud
 very thin, lightly
 curved and rounded tip
KS-1230



Spud & Gouge


In reversible handle
 for removing foreign bodies
KS-1231



Wills Eye Hospital

Foreign body spud & bar
 delicate curved spud on
 one end and small barr
 for rust ring reoval
 on opposite end
KS-1232





West
 Irrigating Cannula
 Blunt tip, 23 gauge
 Side Opening
CS-1233




Kelman
 Cyclodialysis cannula
 spatulated tip
 side opening, 19 gauge
 left opening
 Right opening
CS-1234



Air injection
 Cannula, angled 5mm from tip
 30 gauge
 27 gauge
CS-1235




Sheets
 Cyclodialysis Cannula
 spatulated tip
 one front and two side
 openings, 19 gauge
CS-1236



Viscoelastic
 Aspirating cannula
 Angled 10mm from tip
 22 gauge
 19 gauge
CS-1237



Troutman
 Alpha-Chymotrypsin Cannula
 delicate olive tip
 25 gauge
CS-1238



Rainin
 Air injection cannula
 Spatulated tip
 Angled at 7mm
 30 gauge
 27 gauge
CS-1239



Randolph
 Cyclodialysis cannula
 Curved spatulated tip
 Angled at 13mm
 21 gauge
 19 gauge
CS-1240




Wfelsh
 Flat Olive tip Cannula
 25 gauge
 straight
 left
 right
CS-1241



Aldgerbrush

Corneal Rust
Ring Remover
battery powered
complete
extra burr 0.5mm
extra burr 1.0mm
BB-1242



Castroviejo

Blade Breaker
and holder
with sliding lock
BB-1243



Castroviejo

Blade Breaker and holder
short model
pointed jaws
squared jaws
BB-1244



Swiss Model

Blade Breaker
and holder
spherical tip
BB-1245



Katena


Double Edge Blades
for blade breakers
pk/10 non-sterile
pk/50 sterile
BB-1246



Elliot

Scleral Trephine
1.0mm dia.
1.5mm dia.
2.0mm dia.
Handel
BB-1247






Simcoe

Irrigating-Aspirating Cannula
 Original model, with aspiration
 Through tubing hub and irrigation
 Through luer-lock hub. 23 gauge
 Thin wall, 15mm long curved
 Tubes, with 0.3mm dia. Aspirating
 Port facing into curve and slightly
 Shorter irrigating tip with 45°
 Angled front opening. Designed
 For true Simcoe I/A technique.


IAC-1248



Swiss

Irrigating-Aspirating cannula
 for aspiration through tubing hub
 and irrigation through luer-lock
 hub. 23 gauge thin wall, 15mm long
 curved tubes, 0.3mm dia. Port on left
 side of aspirating tube. Slightly
 angled front opening.


IAC-1249



Reversed simcoe

Irrigating-Aspirating Cannula
 Same as K7-4300 but tubes reversed.
 For surgeons who like the Simcoe I/A
 Cannula, but prefer to use a more
 conventional technique of aspirating
 through the luer-lock hub cannula.


IAC-1250



Ambrosio

Irrigating-Aspirating
 Cannula For aspirating
 through tubing hub and
 Irrigation through luer-lock
 hub. 23 gauge thin wall,
 15mm long, curved,
 front openings.
 smooth tip
 rough tip for polishing

IAC-1251





Kelman
 Cystotome
 With sharp point and
 Sharp cutting edge
 22 gauge
CTS-1252



Kratz
 Cystotome
 With sharp point and
 Sharp cutting edge
 Vaulted, 22 gauge
CTS-1254



Kelman
 Cystotome
 With sharp point and
 Blunt blade, 22 gauge
CTS-1256



Kelman
 Double Blade Cystotome
 With sharp points and
 Blunt blades, 22 gauge
CTS-1258



Double-Edge
 Cystotome
 With sharp point and
 sharp cutting edges
 22 gauge
CTS-1260



Knolle-Kelman
 Cystotome
 small blade
 with sharp point and
 sharp cutting edge
 22 gauge
CTS-1253



McIntyre
 Cystotome
 angled 10mm from tip
 with sharp point, 23 gauge
CTS-1255



Lieppman
 Micro Cystotome
 extra small, with sharp point
 and sharp cutting edge
 22 gauge
CTS-1257



Graether
 Collar Button
 micro iris retractir
 and posterior capsule polisher
 with side irrigation
 straight
 angled
CTS-1259



Shammas
 Posterior Capsule
 Polishing Curette
 1mm diameter cup with
 sandblasted base
 22 gauge
CTS-1261



Lamellar Blade
 sterile, disposable (pk/3)
 #66-60
 #66F
KBH-1263

Miniature Blade
 Sterile, disposable
 #63 (pk/6)
 #64 (pk/12)
 #65 (pk/12)
 #67 (pk/12)
 #69 (pk/6)
 #69B (pk/6)
KBH-1262

Super Blade
 sterile, disposable (pk/12)
 #76-15,15*
 #76-30,30*
 #76-45,45*
KBH-1264

Chuck Handle
 #3H hexagon 10cm
 #3K round 10cm
 #3KA round 7.5cm
 #3KB round 5.0cm
KBH-1265




Stiletto Knife
 For Micro-Vitreo-Retinal
 Stab incisions
 Sterile, disposable (pk/3)
 20 gauge
 19 gauge
KBH-1266



Bard-Parker Handle
 For disposable slot blades
 #3
 #9
KBH-1267




Bard-Parker Blade
 Sterile, disposable (pk/50)
 #11
 #15
KBH-1268





Drews

Posterior Capsule Polisher
 Lightly angled, 2mm I.D loop
 With sandblasted base
 0.2mm irrigating port facing
 into circle, 25 gauge
CP-1269



Kraff

Posterior Capsule
 Polishing Curette
 1.75mm diameter cup
 with semi-sharp edge.
 Irrigating port in center
 Of cup, 25 gauge
CP-1271



Kratz

Scratcher
 Front opening with
 Carbide impregnated tip
 Lightly angled, 22 gauge,
CP-1273



Rubman-katzin

Posterior Capsule Polisher
 Carbide impregnated tip
 With smooth closed end.
 Irrigating side port facing
 Into curve, 22 gauge
CP-1275



Shapiro

Posterior Capsule Polisher
 Carbide impregnated tip
 With dual 0.3mm sideports
 curved. 23 gauge
CP-1277



Jensen

Posterior Capsule Polisher
 carbide impregnated on
 Posterior surface of olive tip
 front opening, gently curved
 25 gauge
 23 gauge
CP-1270



Jensen

Capsule Polisher
 carbide impregnated
 All around olive tip
 front opening, gently curved
 25 gauge
 23 gauge
CP-1272



Knolle

Posterior Capsule Polisher
 2mm angled tip, front opening
 finely sandblasted
 30 gauge
 23 gauge
CP-1274



Simcoe

Cortex Extractor and Capsule Polisher
 side opening facing into curve
 23 gauge, sandblasted
 13mm
 21mm
CP-1276




Hydro-Dissection

Cannula, 25 gauge
 Angled flat tip
CP-1278



Bishope-Harmon
 Irrigating Cannula
 Flat tip, angled at
 8mm from the tip
 sterile, disposable
 19 gauge
CP-1279



Hydrodissection
 Cannula, flat tip
 Angled at 7mm
 From the tip
 Sterile, disposable
 25 gauge
CP-1280



Hydrodissection
 cannula, U-Shaped
 with flat tip
 sterile, disposable
 25 gauge
CP-1282



J-Shaped Cannula
 Sterile, disposable
 Straight, 25 gauge
 angled, 25 gauge
 straight, 23 gauge
 angled, 26 gauge
CP-1284



Simcoe
 Cortex Extractor
 23 gauge X 13mm
 0.3mm aspiration port
 sterile, disposable
CP-1286




Simcoe
 Irrigating-Aspiration
 Cannula, 23 gauge thin wall
 0.3mm aspiration port
 with tubing and adaptor
 sterile, disposable
CP-1281




Kartz
 Capsule Polisher
 end opening, angled at
 7mm from the tip
 sterile, disposable
 23 gauge
 21 gauge
CP-1283



Kartz
 Capsule Polisher
 close tip with side port, angled
 at 7mm from the tip
 sterile, disposable
 23 gauge
 21 gauge
CP-1285



Silicone
 Tip Cannula
 angled at 7mm from
 the tip, 27 gauge
 sterile, disposable
CP-1287





Biploar Forceps

Jeweler's Type
pointed tips, straight
plain
insulated
19 gauge


DI-1288



Biploar Forceps

McPherson Type
lightly blunted tips
straight

DI-1289



Biploar Forceps

McPherson Type
lightly blunted tips
angled

DI-1290



Biploar

Coaptation Forceps
pointed tips


DI-1291



Biploar Cord

disposable(10/box)
reusable

DI-1292





McIntyre-Binkhorst
 Irrigating Cannula
 smooth blunt tip
 front opening, 26 gauge
 right
 straight
 left
C-1293



Gimbel
 Hydro-dissection Cannula
 smooth blunt tip
 front opening, 30 gauge
C-1295



Knolle-Pearce
 Irrigating Vectis
 dual fluid flow through
 three front opening ports
 23 gauge loop
C-1297



Drews
 Irrigating Vectis
 dual fluid flow through
 one front and two backward
 opening ports, 23 gauge loop
C-1299



Sheets
 Irrigating vectis
 dual fluid flow through
 three large front opening
 ports, serrated 20 gauge loop
C-1294




Hofmann
 Nucleus loops
 3mm wide, serrated
 extra delicate
 non-irrigating
C-1296



Simcoe
 Irrigating Lens Nucleus Loop
 two front opening ports
 finely serrated, strong curve
 25 gauge
 right
 left
C-1298



Anis
 Irrigating Vectis
 dual fluid flow through
 two front opening ports
 23 gauges heart shaped loop
 enables surgeon to retract
 iris and capsular flap
 while irrigating
C-1300





Nishi

Nucleus Delivery Cannula
flat insertion plate with
central irrigating port.
Restraining fin limits
depth of insertion and
may be used to apply gentle
pressure to the sclera.
Tip of insertion plate is
placed under the nucleus,
while continuous irrigation
helps to float the nucleus
forward across the plate,
reducing zonular stress.
Right
left
C-1301



Saiz
Irrigating Capsule Retractor
20 gauge, straight
C-1302



Kratz
Lens Manipulating Cannula
0.2mm diameter tip
front opening, 30 gauge
C-1303



Thurmond

**Nucleus Irrigating Cannula
and Iris Retractor, front opening**
smooth downward curved tip
19 gauge
C-1304



Osher

Lens Manipulating Cannula
with sideport and pigtail tip
24 gauge
C-1305



Nishi


Irrigating Capsule Retractors
for hydro-dissection of the
lens nucleus
C-1306



Kuglen


Irrigating Lens Manipulator
disc shaped tip with 0.25mm
diameter peg protruding
from its center, 25 gauge
straight
angled
C-1307






PLANNED EXTRA-CAPSULAR CATARACT EXTRACTION

- | | |
|---|---|
| 1. 2001 Bishop-Harmon A/C Irrigator | 23. 2023 McPherson Tying Forceps, straight |
| 2. 2002 Jones Towel Clamp | 24. 2024 McPherson Tying Forceps, angled |
| 3. 2003 Hartmann Mosquito Forceps, straight | 25. 2025 Kelman-McPherson Forceps |
| 4. 2004 Serrefine, small, straight | 26. 2026 Utrata Capsulorrhexis Forceps |
| 5. 2005 Lieberman Speculum, K-wire | 27. 2027 Sinskey II Lens Hook, angled |
| 6. 2006 Barraquer Wire Speculum, large | 28. 2028 Kuglen Iris Hook, angled |
| 7. 2007 Chuck Handle (2) | 29. 2029 Barraquer Iris Spatula |
| 8. 2008 Super Sharp Blade, 15 deg. (bx) | 30. 2030 Barraquer Needle Holder, curved |
| 9. 2009 Lamellar Blade #66 (bx) | 31. 2031 Castroviejo Needle Holder, straight |
| 10. 2010 Diamond Step Knife | 32. 2032 Lewis Lens Loop |
| 11. 2011 K-Sponge Spear (bx) | 33. 2033 Alfonso Trisector, oval |
| 12. 2012 Castroviejo Corneal Scissors, small, right | 34. 2034 Alfonso Cutting Platform |
| 13. 2013 Castroviejo Corneal Scissors, small, left | 35. 2035 Alfonso Nucleus Forceps |
| 14. 2014 Westcott Tenotomy Scissors | 36. 2036 Alfonso-McIntyre Irrigating Spoon, small |
| 15. 2015 Westcott Stitch Scissors | 37. 2037 Alfonso-McIntyre Irrigating Spoon, large |
| 16. 2016 Vannas Scissors, curved | 38. 2038 Anis Cortex Aspiration Cannula |
| 17. 2017 Stern-Gills Scissors | 39. 2039 Air Injection Cannula, 27ga. |
| 18. 2018 Castroviejo Caliper | 40. 2040 McIntyre A/C Cannula, 26ga. |
| 19. 2019 Bonaccolto Conjunctiva Forceps | 41. 2041 Jensen Polisher, 23ga. |
| 20. 2020 Harms-Colibri Forceps, 0.12mm | 42. 2042 Simcoe I/A Cannula, curved |
| 21. 2021 Bonn Forceps, 0.12mm | 43. 2043 Brierley Chamber Maintainer |
| 22. 2022 Bishop-Harmon Forceps | |



CHANG CONTINGENCY SET


1. 2044 Aspiration Handpiece, 21 gauge
2. 2045 Irrigation Handpiece, 21 gauge
3. 2046 Infusion Cannula, 22 gauge
4. 2047 Corbin Sub-Tenon's Cannula
5. 2048 Simcoe Irrigating Lens Loop, right
6. 2049 Simcoe Irrigating Lens Loop, left
7. 2050 Castroviejo Corneal Scissors, medium, curved



IOL EXPLANTATION SET


1. 2051 Scimeca Haptic Cutter
2. 2052 Stolte IOL Cutter, titanium
3. 2053 Katena IOL Cutter
4. 2054 Sterilization Case, medium, plastic






EXTRA-CAPSULAR CATARACT EXTRACTION USING PHACOEMULSIFICATION

- | | |
|--|---|
| 1. 2055 Bishop-Harmon A/C Irrigator | 22. 2076 Utrata Capsulorrhexis Forceps |
| 2. 2056 Lieberman Speculum, K-wire | 23. 2077 Kelman-McPherson Forceps, angled |
| 3. 2057 Temporal Speculum, open wire | 24. 2078 Ernest-McDonald IOL Inserter |
| 4. 2058 Diamond Knife, tri-facet | 25. 2079 Livernois IOL Folder |
| 5. 2059 Phaco-4 Diamond Knife, 3.0mm, angled | 26. 2080 Katena Needle Holder, curved, w/lock |
| 6. 2060 Chuck Handle | 27. 2081 Ernest Nucleus Cracker |
| 7. 2061 Super Sharp Blade, 15 deg. (bx) | 28. 2082 Iris Retractors, (set 6) in case |
| 8. 2062 Hartmann Mosquito Forceps, str. | 29. 2083 Simcoe I/A Cannula, curved |
| 9. 2063 Jones Towel Clamp | 30. 2084 McIntyre-Binkhorst Cannula, right |
| 10. 2064 Rosen Phaco Splitter | 31. 2085 McIntyre-Binkhorst Cannula, left |
| 11. 2065 Kuglen Iris Hook, angled | 32. 2086 Hydro-Dissection Cannula, 25ga. |
| 12. 2066 Jaffe-Bechert Rotator, 9mm | 33. 2087 Air Injection Cannula, 27ga. |
| 13. 2067 Lester IOL Manipulator, angled | 34. 2088 Jensen Polisher, 23ga. |
| 14. 2068 Westcott Tenotomy Scissors | 35. 2089 Stevens Sub-Tenon's Cannula, 19ga. |
| 15. 2069 Castroviejo Caliper, straight | 36. 2090 Fine-Thornton Fixation Ring, 13mm |
| 16. 2070 Stern-Gills Scissors | 37. 2091 Bipolar Forceps, angled, blunt |
| 17. 2071 Colibri Forceps, 0.12mm teeth | 38. 2092 Bipolar Cord |
| 18. 2072 Castroviejo Forceps, 0.12mm | 39. 2093 McIntyre Reverse Cystotome |
| 19. 2073 Fechtner Conjunctiva Forceps | 40. 2094 Vannas Scissors, curved |
| 20. 2074 Tennant Tying Forceps, straight | |
| 21. 2075 Tennant Tying Forceps, curved | |



PHACO PRECHOP SET

- | | |
|--|--|
| 1. 2095 Lieberman Speculum, V-wire | 12. 2106 Utrata Capsulorrhexis Forceps |
| 2. 2096 Maloney Keratometer, titanium | 13. 2107 Castroviejo Caliper, straight |
| 3. 2097 Ernest Nucleus Cracker | 14. 2108 Westcott Tenotomy Scissors |
| 4. 2098 Sinskey II Lens Manipulator | 15. 2109 McPherson-Westcott Stitch Scissors |
| 5. 2099 Aspiration Handpiece, 21ga. w/ textured tip | 16. 2110 Sauter Hydrodissection Cannula, 27ga. |
| 6. 2100 Irrigation Handpiece, 21ga. dual 0.5mm sideports | 17. 2111 Fukasaku Hydrochop Cannula, 23ga. |
| 7. 2101 Phaco-4 Diamond Knife, 2.8mm, angled | 18. 2112 Hydrodissection Cannula, 25ga. |
| 8. 2102 Diamond Knife, tri-facet, 1mm | 19. 2113 Cystotome, formed, 25ga. (bx/10) |
| 9. 2103 Air Injection Cannula, 30ga. | 20. 2114 Akahoshi Phaco PreChopper |
| 10. 2104 Kratz Scratcher, 22ga. | 21. 2115 Akahoshi Combo PreChopper |
| 11. 2105 Kuglen Hook, angled | |



MICRO INCISION CATARACT SURGERY (MICS) INSTRUMENTS


1. 2116 MICS Diamond Knife, trapezoid, 1.8/2.2mm, angled
2. 2117 Alio-Kelman IOL Hook
3. 2118 Alio MICS Scissors, 23ga.
4. 2119 Alio MICS Capsulorrhexis Forceps, 23 ga.
5. 2120 Giannetti MICS Capsulorrhexis Forceps

6. 2121 Ernest Micro Nucleus Cracker
7. 2122 Irrigating Handpiece, 23ga dual 0.3mm ports
8. 2123 Aspirating Handpiece, 0.3mm w/ textured tip



MANUAL PHACOFRACTURE CATARACT EXTRACTION--AS USED BY EDUARDO ALFONSO, M.D.

- | | |
|---|---|
| 1. 2124 Lieberman Speculum, V-wire | 15. 2138 Bucher-Kansas Nucleus Trisector |
| 2. 2125 Eye Scissors, standard, straight | 16. 2139 Kansas Nucleus Bisector |
| 3. 2126 Westcott Tenotomy Scissors | 17. 2140 Alfonso Nucleus Forceps, 3x4 |
| 4. 2127 Castroviejo Forceps, 0.12mm | 18. 2141 Alfonso-McIntyre Irrigating Spoon, small |
| 5. 2128 Chuck Handle (3) | 19. 2142 Alfonso-McIntyre Irrigating Spoon, large |
| 6. 2129 Lamellar Blade #66 (bx/5) | 20. 2143 McIntyre Infusion Handle |
| 7. 2130 Super Sharp Blade, 15 deg. (bx/5) | 21. 2144 Air Injection Cannula, 27ga. |
| 8. 2131 Phaco Blade 3.0mm, angled (bx/5) | 22. 2145 Anis Cortex Cannula, straight |
| 9. 2132 Castroviejo Caliper, straight | 23. 2146 Brierley Chamber Maintainer |
| 10. 2133 Utrata Capsulorrhexis Forceps | 24. 2147 Bipolar Forceps, straight, pointed |
| 11. 2134 Kuglen Lens Manipulator, angled | 25. 2148 Dodick IOL Folding Forceps |
| 12. 2135 Alfonso Cutting Platform, 2mm | 26. 2149 Buratto Soft IOL Insertor |
| 13. 2136 Alfonso Cutting Platform, 3mm | |
| 14. 2137 Alfonso Nucleus Trisector, oval | |



MANUAL SMALL INCISION CATARACT SURGERY (MSICS) SET--AS USED BY GLENN STRAUSS, MD


1. 2150 Barraquer Wire Speculum, large
2. 2151 Wilder Lens Loop, serrated
3. 2152 Sinskey II Lens Hook, angled
4. 2153 Barraquer Needle Holder, curved, titanium
5. 2154 Harms-Colibri Forceps, 0.12mm
6. 2155 Shepard Lens Forceps
7. 2156 Gaskin Fragment Forceps

8. 2157 Barraquer Iris Spatula, 0.5mm
9. 2158 Hartmann Mosquito Forceps, straight
10. 2159 Vannas Scissors, straight
11. 2160 Westcott Stitch Scissors
12. 2161 Strauss-Simcoe I/A Cannula
13. 2162 Utrata Capsulorrhexis Forceps



LIMBAL RELAXING INCISION SET


1. 2163 LRI Step Knife, 3 Steps (0.5, 0.55, 0.6)
2. 2164 Henderson Alignment Marker
3. 2165 Henderson Degree Gauge
4. 2166 Henderson Limbal Relaxing Incision Marker
5. 2167 Chu Limbal Relaxing Incision Marker
6. 2168 Marking Pen, sterile, disposable (pk/10)



TORIC IOL SET


1. 2169 Henderson Degree Gauge
2. 2170 Henderson Alignment Marker
3. 2171 Henderson Toric IOL Marker
4. 2172 Henderson MICS Capsule Polisher
5. 2173 Marking Pen, sterile, disposable (pk/10)






GLAUCOMA SURGERY

- | | |
|--|---|
| 1. 2174 Bishop-Harmon A/C Irrigator | 16. 2189 Castroviejo Caliper, straight |
| 2. 2175 Jones Towel Clamp | 17. 2190 Diamond Step Knife, 10-facet |
| 3. 2176 Hartmann Mosquito Forceps, straight | 18. 2191 Chuck Handle (2) |
| 4. 2177 Schott Eye Speculum, adjustable | 19. 2192 Lamellar Blade #66 (bx) |
| 5. 2178 Barraquer Wire Speculum, large | 20. 2193 Super Sharp Blade, 15 deg. (bx) |
| 6. 2179 Eye Scissors, straight | 21. 2194 Marking Pen, sterile, disposable (pk/10) |
| 7. 2180 Westcott Tenotomy Scissors | 22. 2195 Luntz-Dodick Trabeculectomy Punch |
| 8. 2181 Barraquer Iris Scissors | 23. 2196 K-Sponge Spear, sterile (bx) |
| 9. 2182 Mini Westcott Scissors | 24. 2197 Blumenthal Conjunctiva Dissector |
| 10. 2183 Katena-Vannas Scissors, 7mm, curved | 25. 2198 Barraquer Needle Holder, curved |
| 11. 2184 McPherson-Westcott Stitch Scissors | 26. 2199 Air Injection Cannula, 27ga. |
| 12. 2185 Fechtner Conjunctiva Forceps (2) | 27. 2200 Harms Trabeculotomy Probe, right |
| 13. 2186 Colibri Forceps, 0.12mm teeth | 28. 2201 Harms Trabeculotomy Probe, left |
| 14. 2187 Bishop-Harmon Forceps | |
| 15. 2188 Fechtner Trabeculectomy Marker | |




CORNEAL TRANSPLANT SURGERY

- | | |
|--|---|
| 1. 2202 Bishop-Harmon A/C Irrigator | 16. 2217 Barron Vacuum Punch, 8.5mm (not shown) |
| 2. 2203 Jones Towel Clamp | 17. 2218 Polack Double Corneal Fixation Forceps |
| 3. 2204 Flieringa Rings (set of 8) | 18. 2219 Castroviejo Forceps, 0.12mm |
| 4. 2205 McNeill Global Fixation Ring | 19. 2220 Castroviejo Forceps, 0.3mm |
| 5. 2206 Barraquer Wire Speculum, large | 20. 2221 Harms-Colibri Forceps, 0.12mm |
| 6. 2207 Lieberman Speculum, solid blade | 21. 2222 Tennant Tying Forceps, straight |
| 7. 2208 Hartmann Mosquito Forceps, straight (2) | 22. 2223 Tennant Tying Forceps, curved |
| 8. 2209 Serrefine, small, straight (2) | 23. 2224 Paton Spatula & Spoon |
| 9. 2210 Troutman-Katzin Scissors, small, right | 24. 2225 Katena Needle Holder, curved w/lock |
| 10. 2211 Troutman-Katzin Scissors, small, left | 25. 2226 Castroviejo Needle Holder, curved w/lock |
| 11. 2212 Westcott Tenotomy Scissors | 26. 2227 Air Injection Cannula, 30ga. |
| 12. 2213 Jaffe Stitch Scissors | 27. 2228 McIntyre A/C Cannula, 26ga. |
| 13. 2214 Barron Vacuum Trephine, 7.5mm | |
| 14. 2215 Barron Vacuum Trephine, 8.0mm (not shown) | |
| 15. 2216 Barron Vacuum Punch, 8.0mm | |


**DEEP ANTERIOR LAMELLAR KERATOPLASTY (DALK) SET**

1. 2229 DALK Corneal Dissector
2. 2230 Rosenwasser DALK Trisector
3. 2231 Bonfadini DALK Spatula
4. 2232 Katzin Corneal Transplant Scissors, right
5. 2233 Katzin Corneal Transplant Scissors, left
6. 2234 DALK Air Injection Cannula, 27ga.
7. 2235 Rosenwasser DALK Cannula, 27ga




ENDOTHELIAL TRANSPLANTATION SET

- | | |
|---|--|
| 1. 2236 Corneal Dissector, straight | 11. 2246 Rubenstein Corneal Roller, sterile (pk/3) |
| 2. 2237 Corneal Dissector, vaulted | 12. 2247 Dahan-Blumenthal A/C Maintainer |
| 3. 2238 Reversed Sinskey Hook | 13. 2248 Barron Vacuum Punch, 8.5mm |
| 4. 2239 Endothelial Stripper | 14. 2249 Barron Vacuum Punch, 8.75mm (not shown) |
| 5. 2240 Rosenwasser Irrigating Endothelial Stripper | 15. 2250 Barron Vacuum Punch, 9.0mm (not shown) |
| 6. 2241 Bores Optic Zone Marker, 9.0mm | 16. 2251 Barron Vacuum Punch, 9.5mm (not shown) |
| 7. 2242 Rosenwasser Lamellar Donor Inserting Shovel | 17. 2252 Barron Artificial Anterior Chamber |
| 8. 2243 Rosenwasser Donor Lamella Inserting Forceps | 18. 2253 Barron Stainless Steel Base |
| 9. 2244 Endothelial Glide Spatula | |
| 10. 2245 DSEK Grasping Forceps, 23ga. | |




LASIK SURGERY

- | | |
|--|---|
| 1. 2254 Lieberman Speculum, K-wire | 10. 2263 Buratto LASIK Flap Forceps |
| 2. 2255 Buratto Flap Protector | 11. 2264 Spatulated LASIK Cannula, 26ga. |
| 3. 2256 Mannis-Buratto Flap Protector | 12. 2265 Buratto LASIK Cannula, 3-ports, 23ga. |
| 4. 2257 Hersh LASIK Spatula | 13. 2266 John IntraLasik Cannula, sterile (pk/10) |
| 5. 2258 Fukasaku LASIK Spatula | 14. 2267 K-Sponge II spears, synthetic (bx) |
| 6. 2259 Katzen LASIK Flap Unzipper | 15. 2268 Marking Pen, standard (pk/10) |
| 7. 2260 Thorlakson Z-LASIK(TM) Flap Lifter | 16. 2269 Lieberman Aspirating Speculum (optional) |
| 8. 2261 Gimbel Stabilization Ring | |
| 9. 2262 Chayet Corneal Marker | |




FEMTOSECOND LASIK SET

- | | |
|--|--|
| 1. 2270 Thorlakson Aspirating Speculum | 6. 2275 Dulaney Corneal Marker |
| 2. 2271 Thorlakson Speculum, wire | 7. 2276 Buratto LASIK Flap Forceps |
| 3. 2272 Buratto LASIK Flap Protector | 8. 2277 Buratto LASIK Cannula, 23ga. |
| 4. 2273 Thorlakson Z-LASIK(TM) Flap Lifter, short handle | 9. 2278 K-Sponge II Spears, synthetic (bx) |
| 5. 2274 IntraLase(R) Flap Lifter, D-E | |

**LASEK SURGERY**


1. 2279 Lieberman Aspirating Speculum
2. 2280 K-Sponge II spears, synthetic (bx)
3. 2281 Alcohol Cannula
4. 2282 Trephine & Well 8/8.5mm

5. 2283 Trephine & Well 9/9.5mm
6. 2284 Micro-Hoe & Epi-Peeler
7. 2285 Flap Repositor




LID SURGERY

- | | |
|--|--|
| 1. 2286 Bishop-Harmon A/C Irrigator | 19. 2304 Desmarres Lid Retractor #3 |
| 2. 2287 Jones Towel Clamp | 20. 2305 Graefe Muscle Hook #2 |
| 3. 2288 Halsted Mosquito Forceps, curved (2) (not shown) | 21. 2306 Fixation Hook Double Sharp, small |
| 4. 2289 Halsted Mosquito Forceps, straight | 22. 2307 Arita Meibomian Gland Compressor |
| 5. 2290 Stevens Scissors, straight | 23. 2308 Bonaccolto Conjunctiva Forceps |
| 6. 2291 Eye Scissors, straight | 24. 2309 Bishop-Harmon Forceps |
| 7. 2292 Chuck Handle | 25. 2310 Castroviejo Forceps, 0.3mm |
| 8. 2293 Performance Miniature Blade #67 | 26. 2311 Helveston Tying Forceps |
| 9. 2294 Westcott Tenotomy Scissors | 27. 2312 Downes Lid Clamp |
| 10. 2295 Westcott Stitch Scissors | 28. 2313 Bipolar Forceps, angled, blunt |
| 11. 2296 Berke Ptosis Forceps | 29. 2314 Bipolar Cord |
| 12. 2297 Ptosis Clamp, 6-pin, concave | 30. 2315 Scleral Shield, small, pair |
| 13. 2298 Barraquer Needle Holder, curved | 31. 2316 Scleral Shield, large, pair |
| 14. 2299 Jaeger Lid Plate, stainless steel | 32. 2317 Helveston Retractor, 11mm |
| 15. 2300 Castroviejo Needle Holder, curved | 33. 2318 Knapp Retractor, 4-prong, sharp |
| 16. 2301 Castroviejo Caliper, straight | 34. 2319 Marking Pen, standard (pk/10) |
| 17. 2302 Desmarres Lid Retractor #1 | |
| 18. 2303 Desmarres Lid Retractor #2 | |



OCULOPLASTIC SET


- | | |
|---|--|
| 1. 2320 Lester-Burch Eye Speculum | 16. 2335 Ptosis Lid Clamp, 6-pin, concave |
| 2. 2321 Schepens Orbital Retractor | 17. 2336 Downes Lid Clamp |
| 3. 2322 Desmarres Lid Retractor #2 | 18. 2337 Francis Chalazion Forceps |
| 4. 2323 Scleral Shield, large (pair) | 19. 2338 Hartmann Mosquito Forceps, straight (2) |
| 5. 2324 Jaeger Lid Plate, stainless steel | 20. 2339 Castroviejo Needle Holder, straight w/ lock |
| 6. 2325 Wright Fascia Needle | 21. 2340 Kalt Needle Holder |
| 7. 2326 Knapp Retractor, 4-prong, blunt | 22. 2341 Freer Periosteal Elevator |
| 8. 2327 Meyerhoefer Curette #3 | 23. 2342 West Lacrimal Chisel, straight |
| 9. 2328 Bunge Evisceration Spoon | 24. 2343 Wilder Lacrimal Dilator #2 |
| 10. 2329 Fixation Hook, double, large | 25. 2344 Bowman Lacrimal Probe #3-4 |
| 11. 2330 Graefe Muscle Hook | 26. 2345 Pigtail Probe |
| 12. 2331 Stainless Steel Ruler | 27. 2346 Bard-Parker Handle #3 |
| 13. 2332 Stevens Scissors, straight | 28. 2347 Bard-Parker Blade #11 |
| 14. 2333 Enucleation Scissors, curved | 29. 2348 Bard-Parker Blade #15 |
| 15. 2334 Tissue Forceps, 1x2, straight | |



PHACO PRECHOP SET

1. 2349 Jaeger Lid Plate, stainless steel
2. 2350 Meyerhoefer Curette #1
3. 2351 Meyerhoefer Curette #2
4. 2352 Meyerhoefer Curette #3
5. 2353 Meyerhoefer Curette #4
6. 2354 Eye Scissors, curved
7. 2355 Manhattan Eye & Ear Suture Forceps


8. 2356 Hunt Chalazion Forceps
9. 2357 Francis Chalazion Forceps
10. 2358 Castroviejo Needle Holder, curved
11. 2359 Miniature Blade #67 (bx)
12. 2360 Chuck Handle



STRABISMUS SURGERY--AS USED BY EUGENE M. HELVESTON, MD


1. 2361 Lieberman Speculum, solid blade, small
2. 2362 Lieberman Speculum, solid blade, large
3. 2363 Hartmann Mosquito Forceps, curved (2)
4. 2364 Serrefine, small, straight (2)
5. 2365 Helveston Retractor, 7mm wide
6. 2366 Helveston Retractor, 9mm wide
7. 2367 Helveston Retractor, 11mm, concave
8. 2368 Helveston Teaser Hook, 6mm (3)
9. 2369 Helveston Finder Hook, large
10. 2370 Helveston Finder Hook, small
11. 2371 Helveston Muscle Hook, 8mm
12. 2372 Helveston Muscle Hook, 10mm
13. 2373 Helveston Muscle Hook, 12mm
14. 2374 Helveston Scleral Ruler

15. 2375 Castroviejo Caliper, straight
16. 2376 Westcott Tenotomy Scissors
17. 2377 Westcott Stitch Scissors
18. 2378 Bonn Forceps, 0.12mm
19. 2379 Moody Fixation Forceps, left
20. 2380 Moody Fixation Forceps, right
21. 2381 Lester Fixation Forceps, 2x3 (3)
22. 2382 Helveston Tying Forceps (2)
23. 2383 Castroviejo Needle Holder, straight, w/lock
24. 2384 Barraquer Needle Holder, curved, w/lock
25. 2385 Bishop-Harmon A/C Irrigator
26. 2386 Chuck Handle
27. 2387 Miniature Blade #69 (bx)




LACRIMAL SURGERY

- | | |
|--|---|
| 1. 2388 Bishop-Harmon A/C Irrigator | 20. 2407 Pigtail Probe w/Suture Holes |
| 2. 2389 Backhaus Towel Clamp (2) | 21. 2408 Crawford Retrieving Hook |
| 3. 2390 Halsted Mosquito Forceps, straight (2) | 22. 2409 Crawford Lacrimal Intubation Set |
| 4. 2391 Halsted Mosquito Forceps, curved (2) (not shown) | 23. 2410 Bowman Lacrimal Probe Set |
| 5. 2392 Kerrison Rongeur, size 1 | 24. 2411 Arita Meibomian Gland Compressor |
| 6. 2393 Nasal Speculum, adult | 25. 2412 Gold Punctal Plug Forceps |
| 7. 2394 Beyer Rongeur | 26. 2413 Dressing Forceps, delicate, straight |
| 8. 2395 Mallet, lead filled stainless steel | 27. 2414 Sauer Fixation Forceps, 1x2 |
| 9. 2396 Performance Miniature Blade #65 | 28. 2415 Bonaccolto Conjunctiva Forceps |
| 10. 2397 Chuck Handle | 29. 2416 Castroviejo Forceps, 0.3mm |
| 11. 2398 Stevenson Retractor 3x3 | 30. 2417 Knapp Retractor, 4-prong, blunt (2) |
| 12. 2399 Shahinian Lacrimal Cannula | 31. 2418 Stevens Scissors, curved |
| 13. 2400 Burnstine Lacrimal Cannula | 32. 2419 Eye Scissors, straight |
| 14. 2401 Lacrimal Cannula, straight | 33. 2420 Castroviejo Needle Holder, curved w/lock |
| 15. 2402 Lacrimal Cannula, curved | 34. 2421 Westcott Stitch Scissors |
| 16. 2403 Castroviejo Lacrimal Dilator | 35. 2422 Tibolt Punctal Scissors |
| 17. 2404 Freer Periosteal Elevator | |
| 18. 2405 West Lacrimal Chisel, straight | |
| 19. 2406 Freer Lacrimal Chisel, curved | |



VITREORETINAL SURGERY


- | | |
|---|---|
| 1. 2423 Fishkind-Castroviejo Speculum | 22. 2444 Castroviejo Forceps, 0.12mm |
| 2. 2424 Serrefine, small, straight (4) | 23. 2445 Castroviejo Forceps, 0.3mm |
| 3. 2425 Halsted Mosquito Forceps, straight (2) | 24. 2446 Bishop-Harmon Forceps (2) |
| 4. 2426 Eye Scissors, straight | 25. 2447 Castroviejo Needle Holder, curved w/lock |
| 5. 2427 Bishop-Harmon A/C Irrigator | 26. 2448 Castroviejo Needle Holder, straight w/lock |
| 6. 2428 Stevens Scissors, straight | 27. 2449 McPherson Tying Forceps, straight |
| 7. 2429 Stevens Scissors, curved | 28. 2450 Helveston Barbie Retractor |
| 8. 2430 Westcott Tenotomy Scissors | 29. 2451 Helveston Retractor, concave, 11mm |
| 9. 2431 Westcott Stitch Scissors | 30. 2452 Viscoelastic Aspiration Cannula, 22ga. |
| 10. 2432 Performance Miniature Blade #64 (bx) | 31. 2453 Randolph Cyclodialysis Cannula, 19ga. |
| 11. 2433 Chuck Handle | 32. 2454 K-Sponge Spears (bx) |
| 12. 2434 Performance MVR Knife, straight, 19ga. | 33. 2455 Bipolar Cord |
| 13. 2435 Schepens Orbital Retractor | 34. 2456 Bipolar Forceps, angled blunt |
| 14. 2436 Graefe Muscle Hook #2 | 35. 2457 Braunstein Caliper, 3.5/4.0mm |
| 15. 2437 Gass Retinal Detachment Hook | 36. 2458 Squeeze Handle Scissors, horizontal 23ga. |
| 16. 2438 Schocket Scleral Depressor | 37. 2459 Squeeze Handle Forceps, serrated 23ga. |
| 17. 2439 Castroviejo Caliper, straight | 38. 2460 Squeeze Handle Forceps, end-gripping 23ga. |
| 18. 2440 Gass Scleral Marker | 39. 2461 Squeeze Handle Forceps, asymmetrical 23ga |
| 19. 2441 Watzke Sleeve Spreading Forceps | 40. 2462 Squeeze Handle Forceps, ILM 23ga. |
| 20. 2442 Bonaccolto Conjunctiva Forceps | |
| 21. 2443 Nugent Utility Forceps, serrated | |



ENUCLEATION SURGERY

1. 2463 Lester-Burch Eye Speculum
2. 2464 Bunge Evisceration Spoon, large
3. 2465 Bunge Evisceration Spoon, small
4. 2466 Wells Enucleation Spoon
5. 2467 Graefe Muscle Hook #2


6. 2468 Stevens Scissors, straight
7. 2469 Enucleation Scissors, strong curve
8. 2470 Tissue Forceps, 1x2
9. 2471 Halsted Mosquito Forceps, curved



PTERYGIUM SURGERY


1. 2472 Miniature Blade #64 (bx)
2. 2473 Chuck Handle
3. 2474 Lancaster Eye Speculum
4. 2475 McPherson Tying Forceps, straight
5. 2476 Castroviejo Forceps, 0.12mm
6. 2477 Gill Corneal Knife
7. 2478 Castroviejo Needle Holder, curved
8. 2479 Katena Ophthalmic Burr, complete

9. 2480 Diamond Burr, 5.0mm dia. ball
10. 2481 Diamon Burr, 5.0mm dia. disc
11. 2482 Eye Scissors, straight, standard
12. 2483 Eye Scissors, straight, small
13. 2484 Bonaccolto Conjunctiva Forceps
14. 2485 Lester Fixation Forceps



OFFICE SET

- | | |
|--|---|
| 1. 2486 Barraquer Wire Speculum, small | 16. 2501 Dressing Forceps, curved |
| 2. 2487 Barraquer Wire Speculum, large | 17. 2502 Arita Meibomian Gland Compressor |
| 3. 2488 Kratz-Barraquer Wire Speculum, large | 18. 2503 Barraquer Cilia Forceps |
| 4. 2489 Schocket Scleral Depressor | 19. 2504 Gold Punctal Plug Forceps |
| 5. 2490 Desmarres Lid Retractor #0 | 20. 2505 Katena Cilia Forceps |
| 6. 2491 Desmarres Lid Retractor #2 | 21. 2506 Jewelers Forceps #3-C |
| 7. 2492 Golf Club Foreign Body Spud | 22. 2507 Suture Removal Forceps |
| 8. 2493 Spud & Gouge | 23. 2508 Francis Chalazion Forceps |
| 9. 2494 Algerbrush II Corneal Burr | 24. 2509 Castroviejo Lacrimal Dilator, D-E |
| 10. 2495 Meyerhoefer Curette #3 | 25. 2510 Bishop-Harmon Bulb & Adaptor |
| 11. 2496 Kimura Platinum Spatula | 26. 2511 Fox Aluminum Eye Shield, adult |
| 12. 2497 New Orleans Lens Loop | 27. 2512 Bishop-Harmon Cannula, 19ga. (bx/10) |
| 13. 2498 Jameson Muscle Hook, small | 28. 2513 Lacrimal Cannula 19ga. (bx/10) |
| 14. 2499 McPherson-Vannas Scissors, curved | 29. 2514 Lacrimal Cannula, 23ga., Straight |
| 15. 2500 Knapp Strabismus Scissors, straight | |



VETERINARY OPHTHALMIC INSTRUMENT SET


- | | |
|---|---|
| 1. 2515 Barraquer Wire Speculum, sm | 16. 2530 Intraocular Scissors, angled |
| 2. 2516 Barraquer Wire Speculum, lg | 17. 2531 Adson Tissue Forceps |
| 3. 2517 Botvin Iris Forceps, 1x2 | 18. 2532 Castroviejo Needle Holder, curved w/lock |
| 4. 2518 Colibri Forceps, 1x2 | 19. 2533 Castroviejo Needle Holder, heavy, curved |
| 5. 2519 Harms-Colibri Forceps, 0.12mm | 20. 2534 Sinskey Lens Hook, angled |
| 6. 2520 Castroviejo Corneal Scissors, sm, right | 21. 2535 Jaffe-Bechert Rotator, 9mm |
| 7. 2521 Castroviejo Corneal Scissors, sm, left | 22. 2536 Chuck Handle (2) |
| 8. 2522 Westcott Tenotomy Scissors | 23. 2537 Bard-Parker Handle |
| 9. 2523 Utrata Capsulorrhexis Forceps | 24. 2538 Jones Towel Clamp (4) |
| 10. 2524 Bishop-Harmon Forceps | 25. 2539 Halsted Mosquito Forceps, straight (2) |
| 11. 2525 Dressing Forceps, delicate, serrated | 26. 2540 Halsted Mosquito Forceps, curved (2) (not shown) |
| 12. 2526 Tissue Forceps, delicate, 1x2 | 27. 2541 Utility Scissors |
| 13. 2527 Chang Phaco Chopper | 28. 2542 Chang Hydrodissection Cannula |
| 14. 2528 Warren Canine Phaco Chopper | |
| 15. 2529 Intraocular Scissors, straight | |




Sterilizing Tray
SC-01417
 Autoclavable
 inclined Scale



Sterilizing Tray
SC-01418
 With Silicon Mat
 (190mmx102mmx38mm)



Sterilizing Tray
SC-01419
 With Silicon Mat
 (254mmx152mmx38mm)



Sterilizing Tray
SC-01420
 With Two Silicon Mat
 (254mmx152mmx38mm)



Our Range of Best Selling Instruments

Ask For it!

info@metech-international.com

*Our Products Are Guaranteed And Meet The Relevant
Worldwide Standards And Our Products Are
Manufactured Entirely Under In Factory.*



METECH
— INTERNATIONAL —

☎ +92 333 7111745 | +92 300 7111745

🖱 www.metech-international.com

✉ info@metech-international.com | export@metech-international.com

📍 Muslim Street, Mubarik Pura, Sialkot, 51310 - Pakistan.

